


## Technical Notes

### Solar Cell



#### Associated Equipment

The following equipment from djb microtech can be used with the Solar Cell:

- Solar Motor L1-1020.00
- Smiley Green Man LED Board L1-1030.00
- Light Meter D3-1010.00
- Decade Resistance Unit A1-1010.00


The mono-crystalline solar panel has a maximum rating of 4.4V, 90mA. It can be easily energised by a bright lamp or sunlight and can be operated in conjunction with our Solar Motor and Smiley Green Man LED Board.

#### Investigations


Listed below are the outlines of some possible investigations to try with your Solar Cell.

- Investigate how the output voltage varies with the area of cell exposed to the light.
- Use filters to investigate how the output voltage varies with colour/wavelength of light.
- Use a light meter and investigate how the output voltage varies with light intensity.
- Determine the optimum electrical load for maximum power transfer for a known light intensity.
- Determine how the maximum power varies with the light intensity when the electrical load is fixed.


#### Circuit for Power Investigations



## Results for Maximum Power Transfer



Maximum power graph with light level = 5800 lux



Maximum power graph with light level = 2200 lux

Graphs drawn using the ALBA Data Logging software - demo version available free from [www.djb.co.uk](http://www.djb.co.uk)

### A challenge

Hunterston B nuclear power station has a maximum power output of 1288MW. If it was to be replaced by the solar cell that you have been using, how many solar cells would be required and what is the area that they would occupy? You can assume the following:

the average daily light intensity is 6,000 lux - averaged over 24 hours.

the cells are packed side by side with no space between them.

To get started you will need to take a lamp and adjust its height so that the light intensity at bench level is 6,000 lux. Now find the maximum power output.

The Technical Notes are available as a coloured pdf in the Teachers section of our website.