## **Physics Workshop**

## **Displacement/Time graph for freely falling object**

- 1. Follow the instructions in the Technical Notes and set up TSA for measuring 'g'
- 2. Release the ball from different heights and note the corresponding time to fall.
- 3. Draw a graph of displacement against time<sup>2</sup>
- 4. Measure the gradient of the graph this should be g/2

## **Additional Experiment**

- In this experiment you will drop a mask through two Light Bridges and calculate the acceleration due to gravity.
- 2. Set TSA to the Acceleration Data mode.
- 3. Connect a light Bridge to each input and align them so that they are one above the other and at least 10cm apart.
- 4. The length of the mask is 10cm. Drop it so that it cuts both IR beams.
- 5. From the data displayed calculate 'g'
- **6.** Once you have done step 3 a few times set TSA to the Acceleration mode and repeat the experiment to show the difference between the two modes.
- Since microcontrollers are very good at timing why do we not get really accurate answers every time? List five sources of error.

## **Apparatus Supplied**

- TSA & plugtop power supply Release Mechanism with ball bearing Timing Plate Leads Mask Metre stick (S)
- 2 Retort stands and boss heads (S)