Physics Workshop

Wavelength of Light

- 1. Connect the plugtop supply to the RayZer Lightbox.
- 2. Play with the lenses, mirrors and prisms.
- 3. Use the diffraction grating (100 l/mm) and project a fringe pattern onto a paper screen placed at the end of the baseboard.
- 4. Calculate the wavelength using $\lambda = dD/L$ where d=slit separation, D=fringe separation, L= distance of screen to grating.
- 5. Repeat with different slit separation.

Apparatus Supplied

RayZer Lightbox + plugtop supply

Lens, Prism & Mirror set

Diffraction grating

Metre Stick (S)

Large white metal screen