djb microtech ltd

Technical Notes Windmill

Associated Equipment

The following djb microtech equipment can be used with the Windmill:

- Solar Motor L1-1020.00
- Smiley Green Man L1-1030.00
- 10F Super Capacitor



The Windmill is an output device and it produces a voltage in the range 1.5V-4.5V.

Things to try with your Windmill

- Connect the Smiley Green Man to your windmill taking care to connect the positive output from the Windmill to the positive terminal on the Smiley Green Man Board. Using a hair drier as a wind source investigate how the brightness of the LEDs varies as the hair drier moves closer to the Windmill you may have to support the Windmill if the wind is too strong.
- Replace the Smiley Green Man with the Solar Motor and repeat the above.
- Connect the output of the Windmill to a voltmeter and repeat the above.
- Turn the base of the windmill in steps of 10 degrees and investigate how the output voltage varies.

djb microtech ltd

Delfie House, 1 Delfie Drive, Greenock, Scotland, PA16 9EN

Phone/Fax: 01475 786540 Email: info@djb.co.uk Website: www.djb.co.uk

Storing Energy using our 10F Super Capacitor



Step 1

Using the hair drier, charge the super capacitor until the voltage is in the range 1.9V - 2.1V. This may take up to 10 minutes. When the voltage is at a suitable value disconnect the lead to the super capacitor with the wind still blowing.

Step 2

Now connect the two leads from the capacitor to the Smiley Green Man board and the LEDs will stay on for **5** minutes.

A challenge

A major UK DIY chain sells roof mounted windmills which have a maximum output of 1kW when the wind speed is 16m/s or greater. Hunterston B nuclear power station has a maximum power output of 1288MW. If Hunsterton B was to be replaced with roof mounted windmills, how many homes would be required? Give two reasons why your answer is too low.

This Technical Note is available as a coloured pdf in the Teachers section of our website.